# Niti-S SPAXUS™ Stent

**Authorized representative in Europe**  
 Représentant autorisé en Europe  
 Bevollmächtigter in Europa  
 Rappresentante autorizzato per l'Europa  
 Representante autorizado en Europa  
 Avrupa yetkili temsilcisi  
 Autoriseret repræsentant i Europa  
 유럽 대리인

**Consult instructions for use**  
 Consulter les instructions d'utilisation  
 Gebrauchsanleitung beachten  
 Consultare le istruzioni per l'uso  
 Consulte las instrucciones de uso  
 Kullanım talimatlarına bakın  
 Konsultér brugsanvisning

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**Catalogue No.**  
 N° de référence  
 Katalog-Nr.  
 N. di catalogo  
 N° de catálogo  
 Katalog No.  
 Referencenummer  
 모델명

**Attention, consult instructions for use**  
 Attention, consulter les instructions d'utilisation  
 Achtung: Gebrauchsanleitung beachten  
 Attenzione, consultare le istruzioni per l'uso  
 Atención, consulte las instrucciones de uso  
 Dikkat! Kullanım talimatlarına bakın.  
 Vær opmærksom, konsultér brugsanvisning  
 주의, 사용상 지시 참조

**Serial No.**  
 N° de série  
 Serien-Nr.  
 N. Seriale  
 N° de serie  
 Seri No.  
 Serienummer  
 일련번호

**Date of Manufacture**  
 Date de fabrication  
 Datum der Herstellung  
 Data di produzione  
 Fecha de fabricación  
 Üretim Tarihi  
 Fabrikationsdato  
 제조일자

**Do not use if package is damaged**  
 Ne pas utiliser si l'emballage est abîmé  
 Nicht verwenden, wenn die Verpackung beschädigt ist!  
 Non utilizzare la confezione se danneggiata  
 No usar si el paquete está dañado.  
 Ambalaj hasar görmüşse tekrar kullanmayın  
 Må ikke bruges, hvis emballagen er beskadiget  
 포장이 파손된 경우 사용금지

**Temperature limitation**  
 Limites de température  
 Temperaturvorgaben  
 Intervallo termico  
 Límites de temperatura  
 Isı sınırlaması  
 Temperaturgrænser  
 보관 온도

**Manufacturer**  
 Fabricant  
 Hersteller  
 Produttore  
 Fabricante  
 Üretici  
 Fabrikant  
 제조사

**Sterilized using ethylene oxide**  
 Stérilisé à l'aide d'oxyde d'éthylène  
 Mit Ethylenoxid sterilisiert  
 Sterilizzato con ossido di etilene  
 Esterilizado por óxido de etileno  
 Etilen oksit kullanılarak sterilize edilmiştir  
 Steriliseret med ethylenoxid  
 E.O Gas 멸균

**Use by(Expiration Date)**  
 A utiliser avant (date d'expiration)  
 Verwendbar bis (Verfallsdatum)  
 Utilizzare entro (data di scadenza)  
 Utilizar antes de (fecha de vencimiento)  
 Son Kullanma Tarihi  
 Bruges før (udløbsdato)  
 유효기간

**Do not reuse**  
 Ne pas réutiliser  
 Nicht wiederverwenden  
 Non riutilizzare  
 No reutilizar  
 Yeniden kullanmayın  
 Må ikke genbruges  
 재사용 금지

**Do not resterilize**  
 Ne pas stériliser à nouveau  
 Nur einmal sterilisieren!  
 Non risterilizzare  
 No reesterilizar  
 Yeniden sterilize etmeyin  
 Må ikke resteriliseres  
 재멸균 금지

**MR Conditional**  
 RM conditionnelle  
 Bedingte magnetische Resonanzsicherheit  
 RM compatibile  
 RM condicional  
 Koşullu MR  
 MR-betinget  
 MR 조건부 입증

**TaeWoong**  
 MEDICAL

## ***MRI Information***

### **MR Conditional**

Niti-S stent was determined to be MR-conditional.

Non-clinical testing and MRI simulations were performed to evaluate the entire family of the Niti-S SPAXUS™ Stent. Non-clinical testing demonstrated that the entire family of the Niti-S SPAXUS™ Stent is MR Conditional. A patient with this device from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only
- Maximum spatial gradient magnetic field of 2,000-gauss/cm (20-T/m)(extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Niti-S SPAXUS™ Stent is expected to produce a maximum temperature rise of 2.1°C after 15-minutes of continuous scanning (i.e., per pulse sequence)

### ***Artifact Information***

In non-clinical testing, the image artifact caused by the Niti-S SPAXUS™ Stent extends approximately 5-mm from this device when imaged using a gradient echo pulse sequence and 3-Tesla MR system.

# English

## User's Manual

### 1. Description

The Niti-S SPAXUS™ Stent consists of an implantable metallic stent and an introducer system.


The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Niti-S SPAXUS™ Stent

**Figure 1. Stent Model**

The Stent is loaded in the introducer system and upon deployment the body of the stent imparts an outward radial force while both flares impart a horizontal inward force, allowing transmural drainage.

Niti-S SPAXUS™ Stent is Fully Silicone Covered and can be removed; (See Warnings).

**Figure 2. Introducer System (Endoscopic)**

- The endoscopic introducer system has a usable length of 180cm

**Endoscopic** Type is recommended

- When approached endoscopically

### 2. Principle of Operation

The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher towards the hub. Retraction of the outer sheath releases the stent.

### 3. Indication for Use

The Niti-S SPAXUS™ Stent is intended for the drainage of a pancreatic pseudocyst or a gallbladder through a transgastric or transduodenal approach.

### WARRANTY

Taewoong Medical Co., Ltd. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not

limited to merchantability or fitness for a particular purpose, with respect to such instruments.

#### **4. Contraindication**

The Niti-S SPAXUS™ Stent is contraindicated for, but are not limited to:

- Hemodynamic instability
- Severe coagulopathy
- All others than indication for use
- Recapturing a stent during its deployment is contraindicated.

#### **5. Warnings**

- The safety and efficacy of this device for use in the vascular system has not been established.
- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Long-term patency of this stent has not been established and regular observation is recommended.
- The stent cannot be resheathed once deployment has been initiated.
- Do not attempt to recapture/reload a stent once its deployment is advanced.

#### **6. Precautions**

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The Niti-S SPAXUS™ Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S SPAXUS™ Stent is intended for single use only. Do not resterilize and/or reuse the device.

#### **7. Instructions in the Event of Damage**

**WARNING:** Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

#### **8. Potential Complications**

Potential complications associated with the use of Niti-S SPAXUS™ Stent may include, but are not limited to:

##### Procedural Complications

- Bleeding
- Stent misplacement
- Inadequate expansion
- Migration
- Pain
- Perforation

### Post Stent Placement Complications

- Bleeding
- Pain
- Perforation
- Stent dislocation
- Stent misplacement or migration
- Stent occlusion
- Fever
- Pancreatitis
- Abscess formation
- Haemorrhage
- Vomiting
- Pneumoperitoneum
- Intraperitoneal leakage
- Peritonitis
- Hematoma
- Inflammation or Infection
- Fistula
- Ulceration
- Sepsis
- Rupture of intracystic artery
- Stent removal failure

### **9. Equipments Required**

- Percutaneous Placement  
Percutaneous placement has not been established yet.
- Endoscopic Placement
  - 0.035" (0.89mm) guidewire at least 450cm long
  - Echoendoscope or therapeutic duodenoscope with instrument channel of 3.7mm or larger
  - Niti-S SPAXUS™ stent
  - 19G FNA needle
  - (Cystotome or needle knife) and/or Dilator (4~6mm balloon dilation catheter)
  - Forceps and/or Snare

### **10. Procedure**

#### **- (50cm introducer usable length)**

Percutaneous procedure has not been established yet.

#### **- (180cm introducer usable length)**

Endoscopy Ultrasound (EUS) should be performed prior to placement of the Niti-S SPAXUS™ Stent to characterize the extent of the lesion and its morphology.

#### **A. Stent Size Determination**

- a) The stent length and diameter should be determined by the doctor after endoscopic and/or fluoroscopic inspection of the lesion.
- b) The length and diameter of the stent should be chosen in a way that both walls of the created transmural drainage remain tightly together, thus preventing any migration.

#### **B. Stent Deployment Preparation**

#### **- Endoscopic Procedure**


- a) Under endoscopic guidance, insert an endoscope until reaching the lesion. Then introduce a needle through the working channel of the endoscope. Advance the needle until reaching the lesion and puncture it.
- b) After puncture of the lesion, insert a guide wire through the needle and advance it across the lesion. Remove the needle slowly and carefully.
- c) After needle removal, insert a dilating device such as a cystotome (or needle knife) and/or balloon catheter along the guide wire until reaching across the lesion and dilate it.
- d) After dilating, remove carefully a dilating device

- e) Remove the stylet from the distal end of the introducer.
- f) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.

### C. Stent Deployment Procedure


**PRECAUTION:** Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent.

- a) Under the fluoroscope and endoscopic guidance, position the introducer system. Inner X-ray marker ('A' of Figure 3) should pass through the wall of a pseudocyst or gallbladder.


**Figure 3**

- b) Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently slide the Y-connector back along the pusher towards the hub.
- d) Under EUS and fluoroscopic guidance, the distal flare is deployed inside of the target site
  - slowly pull back the Y connector when the outer X-ray ring overlaps with the x-ray marker on the inner sheath
  - check the opening of the distal flare


**Figure 4**

- e) Pull back the entire delivery system till seeing the blue marker of outer sheath under endoscopic view

**Figure 5**

- f) Deploy the proximal flare of the stent under endoscopic guidance while making sure that the stent connects both walls together.

**Figure 6**

**CAUTION** Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage.

#### **D. After Stent Deployment**

- Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- Carefully remove the introducer system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal)
- Balloon dilatation inside the stent can be performed if judged necessary.

#### **11. Perform Routine Post-implant Procedures**

- Assess the good position of the stent and effective drainage. A Stent may require up to 1 to 3 days to expand fully.
- Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- Observe the patient for development of any complications.

#### **12. Instructions for Removal of Niti-S SPAXUS™ Stent (See Warnings)**

Grasp the Stent with forceps and/or collapse the proximal end of the Stent with snare then carefully and gently retrieve the Stent.

#### **Reuse Precaution Statement**

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your representative of Taewoong Medical Co., Ltd. For a single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Storage: Store at room temperature (10~40 °C).**

**Disposal Requirements:** The introducer system of Niti-S SPAXUS™ Stent must be properly sealed and disposed in compliance with the regulation of local or hospital at the end of its use.

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