



Endoscopic Ultrasound Radiofrequency Ablation Electrode

Features

- Feasible and safe ablative treatment for soft tissue with a focused application of controlled energy and minimal procedure time.
- Easy to use, all-in-one system, that reduces procedural steps and time of RFA treatment.
- Evident device visualization under EUS.
- Insulated inner cooling system which eliminates charring of soft tissue.
- Controlled RFA environment with specialized VIVA Combo™ generator.



Ablation volume with EUSRA™ 19 Gauge



Ordering Information

Code	Tip Length (mm)	Flexible Length (mm)	Tip Diameter
19-05E	5	1393 + 20	19G (1.06mm)
19-07E	7		
19-10E	10		

- **EUS-guided radiofrequency ablation for management of pancreatic insulinoma by using a novel needle electrode** by Sundeep Lakhtakia et al [Gastrointest Endosc. 2016 Jan;83(1):234-9]
 - **Initial experience of EUS-guided radiofrequency ablation of unresectable pancreatic cancer** by Tae Jun Song et al [Gastrointest Endosc. 2016 Feb;83(2):440-3]
 - **EUS-guided Radiofrequency Ablation (EUS-RFA) of Solid Pancreatic Neoplasm Using an 18-gauge Needle Electrode: Feasibility, Safety, and Technical Success** by Stefano Francesco Crino et al [Gastrointest Liver Dis. 2018 Mar;27(1):67-72]
 - **Technique, safety, and feasibility of EUS-guided radiofrequency ablation in unresectable pancreatic cancer** by Filippo Scopelliti et al [Surg Endosc. 2018 May 15]
 - **Endoscopic ultrasound-guided radiofrequency ablation for management of benign solid pancreatic tumors** by Choi Jun-Ho et al [Endoscopy, 2018 May 4]
 - **Endoscopic ultrasound-guided radiofrequency ablation for pancreatic neuroendocrine tumors and pancreatic cystic neoplasms: a prospective multicenter study** by Marc Barthet et al [Endoscopy. 2019 Jan 22]
 - **EUSRA™ Endoscopic Ultrasound guided Radiofrequency Ablation electrode**
- Caution:** Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE: The EUSRA RF Electrode is indicated for coagulation and ablation of soft tissue when used in conjunction with compatible radio frequency generator.

CONTRAINDICATION: There is a risk that error may result due to the radiofrequency current when it comes to the patients who get transplanted with the pacemakers and cardioverter/defibrillators what are used to put into the body, Thus, do not use the radiofrequency lesion generator and electrode.

POTENTIAL COMPLICATIONS: Burning by overheating of surgical unit; Dangers from inexperienced operator's using; Side effects or cross infection from reuse; Weakness of liver functions; Delayed bleeding in the operated body parts; Recurrence of cancer; Syndromes after RFA treatment such as (perforation, respiratory depression or arrest, abdominal pain, fever, nausea, vomit, right shoulder joint pain and chest discomfort and headache) might occur.

The Potential complications are not limited to the above, and other complications may occur by using Endoscopy.

See Instructions for use and User's guide for full product information.