

Features

Stent Lengths

Multiple lengths (6 cm to 12 cm).

Delivery System

Preloaded in a 10FR delivery system.

Delivery Lengths

Available in 220 cm delivery catheter lengths.

High Conformability


Flexible nitinol design for natural conformity and secure placement.

Endoscopic Marker

Yellow endoscopic marker designed for accurate placement under direct visualization.

Radiopaque Markers


Radiopaque marker: 4(four) at both ends and 3(three) in the middle.



Ordering Information

Colonic Comfort Stent

| Code | Stent | | Delivery | |
|-------------|---------------|-------------|--------------|---------------------|
| | Diameter (mm) | Length (cm) | Profile (fr) | Usable Length (cm)* |
| CDT2206N-22 | 22 | 6 | 10 | 220 |
| CDT2208N-22 | | 8 | | |
| CDT2210N-22 | | 10 | | |
| CDT2212N-22 | | 12 | | |
| CDT2406N-22 | 24 | 6 | 10 | 220 |
| CDT2408N-22 | | 8 | | |
| CDT2410N-22 | | 10 | | |
| CDT2412N-22 | | 12 | | |




Indication for Use:

The Niti-S Colonic Comfort Stent is indicated for the palliative treatment of colorectal strictures produced by malignant neoplasms and to relieve large bowel obstruction prior to colectomy in patients with malignant structures.

Duodenal Comfort Stent

| Code | Stent | | Delivery | |
|-------------|---------------|-------------|--------------|---------------------|
| | Diameter (mm) | Length (cm) | Profile (fr) | Usable Length (cm)* |
| PDT2206N-22 | 22 | 6 | 10 | 220 |
| PDT2208N-22 | | 8 | | |
| PDT2210N-22 | | 10 | | |
| PDT2212N-22 | | 12 | | |



Indications for Use:

The Niti-S Duodenal Comfort Stent is indicated for the palliative treatment of pyloric or duodenal obstructions caused by malignant neoplasms.