

Features

Multiple Scope Compatibilities

180 cm (therapeutic gastroscope compatible) and 220 cm (colonoscope compatible) delivery catheters.

Exclusive Diameter Sizes


18, 20 and 22mm diameter stents and large flanges (26 & 28mm). The 22mm is the largest diameter through the scope stent available in the US to treat a variety of unique patient needs.

100% Handmade with a Single Woven Nitinol

Designed to deliver optimal radial and axial force, each stent is crafted to conform seamlessly to patient anatomy. This craftsmanship ensures every device meets the highest standards, supporting the best possible outcomes.

Wide Range of Lengths

Available in (6 to 15 cm), allowing you to choose the perfect stent to cover any indication.



Ordering Information

Fully Covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)*
EST1806F	18	6	10.5	180
EST1808F		8		
EST1810F		10		
EST1812F		12		
EST1814F		14		
EST1815F	20	15	10.5	180
EST2006F		6		
EST2008F		8		
EST2010F		10		
EST2012F		12		
EST2014F	22	14	10.5	180
EST2015F		15		
EST2206F		6		
EST2208F		8		
EST2210F		10		
EST2212F	22	12	10.5	180
EST2214F		14		
EST2215F		15		
EST2206F-22	22	6	10.5	220
EST2208F-22		8		
EST2210F-22		10		
EST2212F-22		12		
EST2214F-22		14		
EST2215F-22		15		

Partially Covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)*
EST1806B	18	6	10.5	180
EST1808B		8		
EST1810B		10		
EST1812B		12		
EST1814B		14		
EST1815B	20	15	10.5	180
EST2006B		6		
EST2008B		8		
EST2010B		10		
EST2012B		12		
EST2015B	22	15	10.5	180
EST2206B		6		
EST2208B		8		
EST2210B		10		
EST2212B		12		
EST2214B	22	14	10.5	220
EST2215B		15		
EST2206B-22		6		
EST2208B-22		8		
EST2210B-22		10		
EST2212B-22		12		
EST2214B-22		14		
EST2215B-22		15		

Indications for Use:

The Esophageal TTS Stent is intended for use in esophageal strictures caused by intrinsic and/or extrinsic malignant tumors.