

Features and Benefits

Double-Layer Hydrophilic Coating for Low-Friction Navigation


- Designed to reduce friction for smooth guidewire advancement and seamless device exchange.

Flexible Polyurethane Tip for Added Safety

- Soft polyurethane tip helps minimize perforation risk while maintaining flexibility.

Powerful Pushability for Confident Access

- Enlarged core wire delivers excellent pushability to help maintain access through challenging anatomy.



Ordering Information

Code	UDI-DI	Qty.	Type	Wire Diameter (in)	Distal Tip Length (cm)	Total Length (cm)	Tip Diameter (in)			
GW3-45S	8809744979915	1 ea	Straight	0.035	5.5	450	0.025			
	8809744979991	10 ea								
GW3-26S	8809744979922	1 ea				6.0		260		
	8809909403996	10 ea								
GW3-45A	8809744979939	1 ea	Angled		6.0	450				
	8809909404009	10 ea								
GW3-26A	8809744979946	1 ea				6.0			260	
	8809909404016	10 ea								
GW2-45S	8809744979953	1 ea	Straight	0.025	7.5	450	0.021			
	8809909404023	10 ea								
GW2-26S	8809744979960	1 ea				7.0		260		
	8809909404030	10 ea								
GW2-45A	8809744979977	1 ea	Angled		7.0	450				
	8809909404047	10 ea								
GW2-26A	8809744979984	1 ea				7.0			260	
	8809909404054	10 ea								

Indications for Use:

Designed to provide secure access to the gastrointestinal tract for the purpose of passing or exchanging endoscopic accessories during procedures such as ERCP and EUS.