


TaeWoong Guidewire



Reliable Access and Smooth Device Exchange

Designed to provide secure access to the gastrointestinal tract for the purpose of passing or exchanging endoscopic accessories during procedures such as ERCP and EUS.

TaeWoong Medical USA
4801 Wilshire Blvd. Suite 102, Los Angeles, CA 90010
E: info@taewoongusa.com T: 213-944-0503
www.taewoongusa.com

©2026 TaeWoong Medical USA. All Rights Reserved.



Features and Benefits

Double-Layer Hydrophilic Coating for Low Friction Navigation

- Hydrophilic coating reduces resistance during advancement.
- Supports smooth device delivery and exchange. Double-layer hydrophilic coating helps achieve low friction force.

Flexible Polyurethane Tip for Added Safety

- Soft, flexible tip helps reduce perforation risk. Flexible polyurethane tip is designed to prevent perforation.

Powerful Pushability for Confident Access


- Enlarged core wire delivers strong push performance through challenging anatomy. The guidewire offers powerful push ability due to its enlarged core wire design.

GW3-45S
GW3-26S

GW3-45A
GW3-26A

GW2-45S
GW2-26S

GW2-45A
GW2-26A



Code	UDI-DI	Qty.	Type	Wire Diameter (mm)	Distal Tip Length (mm)	Total Length (mm)	Tip Diameter (mm)
GW3-45S	8809744979915	1 ea	Straight	0.78	55	4500	0.64
	8809744979991	10 ea					
GW3-26S	8809744979922	1 ea					
	8809909403996	10 ea					
GW3-45A	8809744979939	1 ea	Angled	0.78	60	4500	
	8809909404009	10 ea					
GW3-26A	8809744979946	1 ea					
	8809909404016	10 ea					
GW2-45S	8809744979953	1 ea	Straight	0.61	75	4500	0.53
	8809909404023	10 ea					
GW2-26S	8809744979960	1 ea					
	8809909404030	10 ea					
GW2-45A	8809744979977	1 ea	Angled	0.61	70	4500	
	8809909404047	10 ea					
GW2-26A	8809744979984	1 ea					
	8809909404054	10 ea					