




# Biliary Slim™ M Stent

6Fr delivery system with tapered tip



## Features

Radiopaque marker: 3 (three) at both ends & 2 (two) in the middle  
Hook and Cross for conformability

DISTRIBUTED BY:  
TaeWoong Medical USA  
4801 Wilshire Blvd, Unit 280 Los Angeles, CA 90010 USA

Email: [info@taewoongusa.com](mailto:info@taewoongusa.com)  
[www.TaeWoongUSA.com](http://www.TaeWoongUSA.com)





## Ordering Information

### Biliary Slim M Stent

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)*
BN0604-6	6	4	6	180
BN0605-6		5		
BN0606-6		6		
BN0607-6		7		
BN0608-6		8		
BN0609-6		9		
BN0610-6		10		
BN0612-6		12		
BN0804-6	8	4	6	180
BN0805-6		5		
BN0806-6		6		
BN0807-6		7		
BN0808-6		8		
BN0809-6		9		
BN0810-6		10		
BN0812-6		12		

#### Indications for Use:

The Niti-S Biliary Slim M Stent is indicated for the palliation of malignant strictures in the biliary tree.