



EUS Biopsy Needles for Enhanced Diagnostic Accuracy



Enhanced Visibility with Spiral Echo Pattern

- Spiral Echo design provides superior visibility under ultrasound guidance
- Spiral grooves along the needle surface minimize reflection dispersion and deliver clearer echogenicity



Spiral Grooves

Needle Options for Enhanced Performance

FNA(Fine Needle Aspiration)

- Back cut bevel design at the distal tip reduces resistance and enhances tissue penetration

FRANSEEN

- Crown-shaped needle design ensures higher tissue yield and better diagnostic reliability

FRANSEEN

FNA (A type)



Spiral Echo Pattern

Ordering Information

Needle Type	Code	Gauge	Needle Material	Working Length
FNA Type (Aspiration)	FM-CTA011	19G	Stainless Steel	145cm
	FM-CTA012	22G		
	FM-CTA013	25G		
	FM-CTA011N	19G	Nitinol	
	FM-CTA012N	22G		
	FM-CTA013N	25G		
Fransen Type	FM-CTA002F	22G	Stainless Steel	

Indications for Use:

The ClearTip FNA and FNB Types are intended for Ultrasonically Guided Fine Needle Aspiration, (FNA) of submucosal and extraluminal lesions of the Gastrointestinal Tract, (e.g., lymph nodes, abnormal tissue in the mediastinum).