

# TaeWoong Cystotome





TaeWoong Medical USA  
4801 Wilshire Blvd. Suite 102, Los Angeles, CA 90010  
E: [info@taewoongusa.com](mailto:info@taewoongusa.com) T: 213-944-0503  
[www.taewoongusa.com](http://www.taewoongusa.com)

©2026 TaeWoong Medical USA. All Rights Reserved.



# TaeWoong Cystotome



## Ordering Information

Model	Electrode Distal Tip OD	Sheath OD	Usable Length	Recommended G/W
OCT1906	1.4mm	2.0mm	1,850mm	0.025" ~ 0.035"
OCT1907	1.4mm	2.3mm		
OCT1908	1.7mm	2.8mm		

### Indications for Use:

The Cystotome is intended for use as an electrosurgical accessory to cannulate the transgastric or transduodenal wall and into a pancreatic pseudocyst that is visibly bulging into the gastrointestinal tract. This device is supplied sterile and is intended for single use only.