


Enhanced precision for complex access

Cystotome is designed to deliver reliable, precise access during EUS-guided procedures — when it matters most. Engineered to advance over a guidewire after the initial FNA puncture, enabling a smooth transition from diagnostic to therapeutic access.



Ordering Information

Code	Electrode Distal Tip OD	Sheath OD	Usable Length	Recommended C/W
OCT1906	1.4mm	2.0mm	1,850mm	0.025" ~ 0.035"
OCT1907	1.4mm	2.3mm		
OCT1908	1.7mm	2.8mm		

Diverse tip size options

6 Fr

OD 2.0mm

7 Fr

OD 2.3mm

8.5 Fr

OD 2.8mm

Indications for Use: The Cystotome is intended for use as an electrosurgical accessory to cannulate the transgastric or transduodenal wall and into a pancreatic pseudocyst that is visibly bulging into the gastrointestinal tract. This device is supplied sterile and is intended for single use only.