

Optimos

# Cystotome





DISTRIBUTED BY:  
TaeWoong Medical USA  
4801 Wilshire Blvd, Suite 280  
Los Angeles, CA 90010 USA

Email: [info@taewoongusa.com](mailto:info@taewoongusa.com)  
[www.TaeWoongUSA.com](http://www.TaeWoongUSA.com)



# Optimos Cystotome



Ordering Information				
Models	1	2	3	Recommended Guide wire
	Sheath O.D.	Tip O.D.	Usable length	
OCT1906	Ø2.0	Ø1.4	185 cm	0.025~0.035"
OCT1907	Ø2.4	Ø1.7		
OCT1908	Ø2.8	Ø1.7		

**Indications for Use:**

The Optimos Cystotome is intended for use as an electrosurgical accessory to cannulate the transgastric or transduodenal wall and into a pancreatic pseudocyst that is visibly bulging into the gastrointestinal tract. This device is supplied sterile and is intended for single use only.