




# Colonic and Duodenal Comfort Stents




TaeWoong Medical USA  
4801 Wilshire Blvd. Suite 280, Los Angeles, CA 90010  
E: [info@taewoongusa.com](mailto:info@taewoongusa.com) T: 213-944-0503  
[www.taewoongusa.com](http://www.taewoongusa.com)

©2025 TaeWoong Medical USA. All Rights Reserved.



# Niti-S Colonic and Duodenal Comfort Stents



## Flexible nitinol design for natural conformity and secure placement

- Designed for conformability to the natural curves of the anatomy.
- Yellow endoscopic marker designed for accurate placement under direct visualization.
- Radiopaque marker: 4(four) at both ends and 3(three) in the middle.

## Ordering Information

### Colonic Comfort Stent

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)*
CDT2206N-22	22	6	10	220
CDT2208N-22		8		
CDT2210N-22		10		
CDT2212N-22		12		
CDT2406N-22	24	6	10	220
CDT2408N-22		8		
CDT2410N-22		10		
CDT2412N-22		12		

#### Indication for Use:

The Niti-S Colonic Comfort Stent is indicated for the palliative treatment of colorectal strictures produced by malignant neoplasms and to relieve large bowel obstruction prior to colectomy in patients with malignant structures.

### Duodenal Comfort Stent

PDT2206N-22	22	6	10	220
PDT2208N-22		8		
PDT2210N-22		10		
PDT2212N-22		12		

#### Indications for Use:

The Niti-S Duodenal Comfort Stent is indicated for the palliative treatment of pyloric or duodenal obstructions caused by malignant neoplasms.